


## Understanding Life One Protein at a Time ...



Notes

It is recommended that the product is reconstituted with sterile water into a final concentration of 1 mg/ml. Store the reconstituted product in aliquots at -20°C. Avoid multiple freeze-thaw cycles and exposure to frequent changes in temperature.

The use of strong acids and bases, strong oxidants, and high concentrations of organic solvents should be avoided to prevent denaturation.

For research use only!



Tel: 1-732-444-1618 • Fax: 1-888-791-1618 • Email: info@lifetein.com